

Weighing room


The Mettler Toledo XP2U Ultra Micro Balance is used to weigh very small samples, even a single microfossil, with optimum measuring performance (stabilization time) in seconds. The microbalance has a capacity to 6g and readability down to 0.1 μg



It poses a fully automatic time- and temperature-controlled internal adjustment and linearization. It has color touchscreen and infrared sensors for hands-free operation. The LabX software connects the microbalance operating system to the PC.



The microbalance is also equipped in an anti-static kit that removes electrostatic charging which has adverse effects on weighing.



Contact person for lab:

- [Tine Lander Rasmussen](#)